



Canine Tests – Price List

Viral diseases

Code	Test name	Specimen	Reporting time	Price LKR
CPV	Canine Parvovirus	Two (2) rectal Swabs	36 Working Hours	3000
CDV	Canine Distemper	1-2 ml of Blood for serum / Blood for plasma/CSF or Swabs	48 Working Hours	3500

Bacterial and other diseases

Code	Test name	Specimen	Reporting time	Price LKR
CB	Canine Babesiosis	Blood for Plasma (1-2 ml)	36 Working Hours	3000
ELS	Canine Ehrlichiosis	Blood for Plasma (1-2 ml)	36 Working Hours	3500

Cattle Test – Price Information

DNA Typing

Code	Test name	Specimen	Reporting time	Price LKR
CaPTT	Cattle Parentage Testing	1-2ml of whole blood	2 weeks	10000

Note:

- 'Sample Collection Guide' on the next page provides information on sample collection
- Samples can be handed over to Genetech from Monday to Friday 8.30 am to 4.30 pm
- The sample will be analysed once the full payment is made.

For more information, please contact us.

Email: info@genetechsrilanka.com

Phone No.: +94112696992/+94112696881

Address: 54 Kitulwatte Road, Colombo 08, Sri Lanka



Sample Collection Guide

Sample/Specimen	Collection	Transport	Storage
Whole blood	1-2 ml of blood in ACD or EDTA (imp- not in Heparin) delivered within 6 hours of collection	Up to 6 hours at room temperature At +4°C (on ice) if transported for more than 6 hours	At +4°C (on ice) for up to 7 days
Clotted Blood	1-2ml of blood delivered within 6 hours of collection	Up to 6 hours at room temperature At +4°C (on ice) if transported for more than 6 hours	At +4°C (on ice) for up to 7 days
Swabs (Nasal, Vaginal, Throat, etc)	Collected into a sterile container and delivered within 6 hours	Up to 6 hours at room temperature At +4°C (on ice) if transported for more than 6 hours	At +4°C (on ice) for up to 7 days At -70°C for more than 7 days
Stool / Urine	Collected into a sterile container and delivered within 6 hours	Up to 6 hours at room temperature At +4°C (on ice) if transported for more than 6 hours	At +4°C (on ice) for up to 7 days At -70°C for more than 7 days
Body Fluids (CSF, Bronchoalveolar Lavage, etc)	Collected into a sterile container and delivered within 6 hours (If contaminated with blood add ACD or EDTA)	Up to 6 hours at room temperature At +4°C (on ice) if transported for more than 6 hours	At +4°C (on ice) for up to 7 days At -70°C for more than 7 days
Bone Marrow	Collected into a sterile container having ACD or EDTA and delivered within 6 hours	Up to 6 hours at room temperature At +4°C (on ice) if transported for more than 6 hours	At +4°C (on ice) for up to 7 days