



## Canine Tests – Price List

### Viral diseases

Code	Test name	Specimen	Reporting time	Price LKR
CPV	Canine Parvovirus genome detection by PCR	Two (2) rectal Swabs	36 Working Hours	4900
CDV	Canine Distemper virus genome detection by PCR	1-2 ml of Blood for serum / Blood for plasma/CSF or Swabs	48 Working Hours	6800

### Bacterial and other diseases

Code	Test name	Specimen	Reporting time	Price LKR
CB	Canine Babesiosis detected by PCR	Blood for Plasma (1-2 ml)	36 Working Hours	4900
ELS	Canine Ehrlichiosis detected by PCR	Blood for Plasma (1-2 ml)	36 Working Hours	4900

## Cattle Test – Price Information

### DNA Typing

Code	Test name	Specimen	Reporting time	Price LKR
CaPTT	Cattle Parentage Testing	1-2ml of whole blood	2 weeks	15000

## Psittacine Beak and Feather Disease (PBFD)

Code	Test name	Specimen	Reporting time	Price LKR
PBFD	Beak and feather disease virus detection by PCR	Feather	2 days	7000

#### Note:

- Sample Collection Guide on the next page provides information on sample collection
- Samples can be handed over to Genetech from Monday to Friday 8.30 am to 4.30 pm
- The sample will be analysed once the full payment is made.

**For more information, please contact us.**

Email: [info@genetechsrilanka.com](mailto:info@genetechsrilanka.com)

Phone No.: +94112696992/+94112696881

Address: 54 Kitulwatte Road, Colombo 08, Sri Lanka



## Sample Collection Guide – Canine Testing

Sample/Specimen	Collection	Transport	Storage
Whole blood (Blood for plasma)	1-2 ml of blood in ACD or EDTA (imp- not in <b>Heparin</b> ) delivered within 6 hours of collection	Up to 6 hours at room temperature  At +4°C (on ice) if transported for more than 6 hours	At +4°C (on ice) for up to 7 days
Clotted Blood (Blood for serum)	1-2ml of blood delivered within 6 hours of collection	Up to 6 hours at room temperature  At +4°C (on ice) if transported for more than 6 hours	At +4°C (on ice) for up to 7 days
Swabs (Nasal, Rectal, Throat, etc.)	Collected into a sterile container and delivered within 6 hours	Up to 6 hours at room temperature  At +4°C (on ice) if transported for more than 6 hours	At +4°C (on ice) for up to 7 days  At -70°C for more than 7 days
Body Fluids (CSF, etc.)	Collected into a sterile container and delivered within 6 hours (If contaminated with blood add ACD or EDTA)	Up to 6 hours at room temperature  At +4°C (on ice) if transported for more than 6 hours	At +4°C (on ice) for up to 7 days  At -70°C for more than 7 days

## Sample Collection Guide – Cattle Testing

Sample/Specimen	Collection	Transport	Storage
Whole blood (Blood for plasma)	1-2 ml of blood in ACD or EDTA (imp- not in <b>Heparin</b> ) delivered within 6 hours of collection	Up to 6 hours at room temperature  At +4°C (on ice) if transported for more than 6 hours	At +4°C (on ice) for up to 7 days